

# DNA Sequencing techniques





1. The “Sanger” Technique
2. 2<sup>nd</sup> generation sequencing (Illumina <sup>TM</sup>, 454 <sup>TM</sup>, Solid<sup>TM</sup>)
3. 3<sup>rd</sup> generation sequencing (Pacific Bio, Oxford Nanopore)

# “Sanger” DNA sequencing

Based on two components

1. Electrophoresis on **polyacrylamide gel**
2. Use of **dideoxyribonucleotides** (ddNTPs)




# Sanger DNA sequencing: Principles




3 major differences with a « normal » PCR :

1. **One single** primer (oligonucleotide)
2. **Robust ADN polymerase** (ex : Sequenase)
3. **NTP mix** → **dNTPs** (large quantity) and **ddNTPs** (at lower concentration)


# Sanger sequencing: Mechanism




# Radioactive “Sanger” sequencing: Mechanism



# Radioactive “Sanger” DNA sequencing: Mechanism




A T G C • Same reaction is repeated for each nucleotide (ddTTP, ddGTP, ddCTP)



- Each reaction will show the position of the complementary dideoxyribonucleotides along the strand
- The upper part of the gel represents the extremity of ...
  - ? of the template
  - ? of the strand that has been synthesised

# Fluorescent “Sanger” sequencing: Mechanism

- Mechanism = the same
- The ddNTPs are labeled with different fluorophores (NOT the dNTPs!)
- Automation! → DNA “sequencer”



# Alternatives to Sanger sequencing?



## 2<sup>nd</sup> generation DNA sequencing


- Pyrosequencing (*454<sup>TM</sup>*)
- Polony sequencing (PCR colonies) (*Illumina<sup>TM</sup>*)

## 3<sup>rd</sup> generation sequencing (**Not yet in the market**)

mostly based on single molecule sequencing



# Pyrosequencing (454 <sup>TM</sup>)




Based on Total DNA sequencing

Follows 4 major steps


1. Break Total DNA (or cDNA) in small fragments
2. Ligate adapters to each newly produced fragment
3. Amplification of fragments by PCR
4. Parallel sequencing of each newly produced fragment **by luminescence (Pyro)**

# How to prepare a 454™ library




1. Break Total DNA (or cDNA) → Produce small fragments

A. DNA Extraction




B. Nebulization (or Sonication)

Small DNA  
fragments  
**Double –  
stranded**  
~ 400 bp




# Pyrosequencing (454™)




## 2. Adapter's ligation

A. Double stranded fragments




Ligase



Adapters



B. Two adapters are ligated and denatured



# The adapters




DNA fragment




**Attachment on streptavidin beads**

# Emulsion PCR

Fragments with adapters  
+  
Streptavidin beads



One fragment per bead



Mixed in a solution that contains Water, Oil, and all the reagents necessary to perform a PCR reaction

Water will form bubbles → PCR MICROREACTORS


# Emulsion PCR




**Each fragment is amplified  
on each bead (within each  
bubble)**



**Results in the cloning of  
millions of fragments, all  
attached on the same bead!**




# Pyrosequencing





**The beads are added  
into “micro-wells”**

**One bead per well!**

**The different beads (i.e.  
the different  
fragments) will be  
sequenced in parallel**





# Pyrosequencing : process






# Pyrosequencing : process



The DNA polymerase incorporates one complementary nucleotide to the fragments attached to one one bead


→ Release of bi-phosphate



The bi-phosphate is used by the Sulfurylase to produce l'ATP

The Luciferase uses ATP to produce LIGHT (PYRO)



# Pyrosequencing : process



The Apyrase degrades the nucleotides that have not been incorporated

The cycle continues, one nucleotide at a time!

The luminescent signal will be proportional to the number of nucleotides that have been incorporated



# Pyrosequencing requires highly specialized instrumentation


Genome sequencer GS10 (Roche)

Utility of Pyrosequencing :

1. Whole genome “de-novo” sequencing
2. Genome re-sequencing
3. Transcriptome sequencing
4. Parallel genotyping




# PCR colony sequencing (Illumina™)



Based on sequencing  
of **total DNA**,  
fragmented *a priori*

Also requires the  
**ligation of specific  
adapters**


Sequencing done on a  
solid surface



# PCR colony sequencing (Illumina™)

## 6 major steps


1. Total DNA (or cDNA) broken into very small fragments
2. **Double - Ligation of adapters** to each, newly produced fragment
3. Selection of fragments with a specific size (200bp)
4. Fragments attached on a **solid** surface
5. Solid-phase PCR → production of PCR colonies (colonies)
6. Fragment denaturation and sequencing **by fluorescence**



# Illumina™ library production


## 1. Total DNA (or cDNA) broken into small fragments

### A. DNA extraction




### B. Nebulization (or Sonication)

Small **double-stranded** DNA fragments  
~ 200 pb




# Illumina™ library production



## 2. Ligation of “UNIVERSAL” Illumina adapters




3. Ligated fragments of a certain size  
are selected




4. Addition of **2 new adapters** by PCR

# New adapters added by PCR








# Attachment of fragments on solid surface (Glass)



# Illumina™ Sequencing






Clusters amplifies the  
Fluorescent signal



# Illumina™ sequencing

One color for each  
nucleotide  
incorporated in the  
different “clusters”



Requires highly  
specific equipment,  
reagents and software



# Comparing sequencing technologies that are readily available



## “Sanger” sequencing :

- Maximum sequence length = 1000 nucleotides
- One single DNA fragment sequenced at a time
- Requires prior knowledge of the sequence to design primers

## Pyrosequencing (454™):

- Maximum sequence length = 400-600 nucleotides
- Hundreds of thousands of fragments sequenced in parallel
- Useful to sequence large chunks of the genome or transcriptomes
- Doesn't require prior knowledge about the genome sequence, content and structure

## PCR Colony sequencing (Illumina™):

- Maximum sequence length = 100 nucleotides
- Millions of fragments sequenced in parallel
- Useful to sequence large chunks of the genome or transcriptomes
- Doesn't require prior knowledge about the genome sequence, content and structure

# Other available Technologies...

SOLID™ sequencing (Sequencing by Oligonucleotide Ligation and Detection - Applied Biosystems):

- Library production = in-between 454 and Illumina
- Very “special” way of sequencing each fragment
- Unique way of reading the sequencing results.

HELICOS™ sequencing:

- No PCR involved!
- True, single DNA molecule sequencing. Each fragment is directly sequenced in parallel, without the need to clone it

# 3<sup>rd</sup> (4<sup>th</sup>??) generation DNA sequencing



## Pacific Biosciences

- Single molecule DNA sequencing
- One DNA molecule + One polymerase!


## Oxford Nanopore

- Single molecule DNA sequencing
- One DNA molecule + One polymerase!

# Genome Assembly

## Shotgun method

### A. DNA extraction



**ADN  
Total**



### B. DNA breakage



### C. Genomic library

**Plasmidic or chromosomic vector  
or with phage vector**




**Library on beads**



**Library on solid phase**




# The « Shotgun » method and vector library




Complete genome:  
2 Mb (very small)



DNA broken into  
1.0-2.0 kb fragments



Insertion into  
plasmids




Random sequencing of clones  
with universal primers


Library including  
overlapping clones




# Sequence assembly



- The different sequences obtained have to be assembled into larger fragments
- We use the « overlap » in sequence between different clones




# Genome assembly – the contigs




- The assembly of sequence fragments into contigs is done using **bioinformatics** (*in silico*)
- Contigs = same concept for all sequencing technologies

# Genome assembly – the contigs



**Contig 1 = 13Kb**  
= region X of the genome




**Contig 2 = 50 kb**  
= region Y of the genome

**Why the contigs dont include entire chromosomes?**

1. Some regions are absent from our library – they have not been included in the vector (probability)
2. Some regions are highly repeated (99% of the cases) = PUBLIC ENEMY NUMBER ONE (...for assemblers!)

# Genome assembly and repeated regions



Repeated regions create identical contigs!




Software for assemblies don't know where to place them

Many regions are repeated in TANDEM → even worse than normal repetitions!

# Genome assembly and repeated regions


1. PCR → with all possible combinations of primers




# Genome assembly and repeated regions

2. Clone and sequence much longer regions of the genome  
→ Increase the probability that a repeated region is already included in the clone

→ Increase the chances of having an overlap



Contig 1 = 13kb  
= region X of the genome




Contig 2 = 50 kb  
= region Y of the genome

.....→ = 2 kpbs = not long enough to link the contigs *in silico*

# Genome assembly and repeated regions

2. Clone and sequence much longer regions of the genome  
→ Increase the probability that a repeated region is already included in the clone

→ Increase the chances of having an overlap



Contig 1 = 13kb

= region X of the genome

Contig 2 = 50 kb

= region Y of the genome

.....➤ = 6 kpbs = long enough to link the contigs *in silico*

# Genome assembly and short reads (454<sup>TM</sup>, Illumina<sup>TM</sup>)

## 3. Paired-Ends or Mate-Pairs:

→ Gives information on the physical distance between **2 sequences** of 2<sup>nd</sup> génération (454 ou Illumina)




The different fragments are “tagged” and sequenced from both extremities.

→ Distance between the sequences obtained is known, because we select fragments of a given size before sequencing (ex: 200pbs)




# Closing a sequence gap




Fragments from different types of libraries  
(plasmid, bacteriophage, 454, Illumina, etc....)


# Genome assembly and short reads (454<sup>TM</sup>, Illumina<sup>TM</sup>)





↓ We know these sequences are separated by X base paires ↓



← → Represents an unknown sequence (NNNN), that we know is located at this position in the genome



# 2000 PCRs later... (+ analyses *in silico*)



# Genomic Survey of the Non-Cultivable Opportunistic Human Pathogen, *Enterocytozoon bieneusi*

Donna E. Akiyoshi<sup>1\*</sup>, Hilary G. Morrison<sup>2\*</sup>, Shi Lei<sup>1</sup>, Xiaochuan Feng<sup>1</sup>, Quanshun Zhang<sup>1</sup>, Nicolas Corradi<sup>5</sup>, Harriet Mayanja<sup>3</sup>, James K. Tumwine<sup>4</sup>, Patrick J. Keeling<sup>5</sup>, Louis M. Weiss<sup>6</sup>, Saul Tzipori<sup>1</sup>

Research

Open Access

**Draft genome sequence of the *Daphnia* pathogen *Octosporea bayeri*: insights into the gene content of a large microsporidian genome and a model for host-parasite interactions**

Nicolas Corradi<sup>\*</sup>, Karen L Haag<sup>†‡</sup>, Jean-François Pombert<sup>\*</sup>, Dieter Ebert<sup>‡†</sup> and Patrick J Keeling<sup>‡\*</sup>

ARTICLE

Received 1 Jun 2010 | Accepted 25 Aug 2010 | Published xx xxx 2010

DOI: 10.1038/ncomms1082

The complete sequence of the smallest known nuclear genome from the microsporidian *Encephalitozoon intestinalis*

Nicolas Corradi<sup>1\*,†</sup>, Jean-François Pombert<sup>1\*,</sup>, Laurent Farinelli<sup>2</sup>, Elizabeth S. Didier<sup>3</sup> & Patrick J. Keeling<sup>1</sup>

Complexity of the assembly